

## Eye Shield

Eye shield is used for a short duration to protect the patient’s eye from accidental injury or rubbing after a cataract surgery.



### Product Description

Product Code	Description	Qty. per box
MIPL/O3/01	Clear polycarbonate Eye shield	100 Pcs. (Unsterile)
MIPL/O3/02	Polypropylene Eye shield	100 Pcs. (Unsterile)

### Features

- Made from high quality clear polycarbonate or polypropylene
- The eye shield covers the patient's eye using an adhesive tape or a patch
- Smooth edges for maximum comfort
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Eye Shields are supplied both non sterile and sterile as per buyers requirement



**MADHU INSTRUMENTS PVT**

A-260, Okhla Industrial Area, Phase-1, New Delhi -110020, India

T: (+91-11) 4270 1030, 31 | E: Sales@madhuinstruments.com | W: www.madhuinstruments.com

Customer Care No. +91 9810866176