



Canabrava Gupta Segment

Canabrava Gupta Segment (Capsule Support Segment) is a partial ring of clear PMMA covering approximately one quadrant with a specially designed scleral fixation arm for permanent fixation.



DR. SERGIO CANABRAVA
BRAZIL
sergiocanabrava@hotmail.com



DR. SATISH C. GUPTA
INDIA
scgupta.dr@gmail.com

Product Description

Product Code	Description					Qty. per box
	Model	Design	Size	Eyelets	Material	
MIPL/D11/CGS	Capsule Support Segment (Canabrava Gupta Segment)		6.00 mm	3	Clear PMMA	1 pc. (Sterile)

Advantages

- Blunt-tipped eyelets to aid in the insertion of the ring.
- Additional eyelet to suture the ring to scleral tissue for better stability of the capsular bag.
- Special angle of the scleral fixation arm avoids tilting during the surgery.
- The eyelet is projecting anteriorly which allows it to sit outside the capsular bag.
- Due to the small size of the segment it is easier to insert inside the bag.

Indications of Canabrava Gupta Segment

- Cases with extreme zonular weakness.
- Cases with missing zonules larger than 4 clock hours
- Lens subluxation

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)