


## Capsular Tension Ring

Capsular tension ring allow expansion and stabilization of the capsular bag to facilitate a safe IOL implantation and maintain good IOL centration. Capsular tension ring is a highly flexible clear PMMA ring of haptic thickness 0.20 mm.



### Product Description

Product Code	Model	Design	Normal Size	Compressed Size	Eyelets	Material	Qty. per box
MIPL/D2/T3	CTR Type 3		12.0 mm	10.0 mm	2	PMMA	1 pc. (sterile)
MIPL/D2/T4	CTR Type 4		13.0 mm	11.0 mm	2	PMMA	1 pc. (sterile)
MIPL/D2/T5	CTR Type 5		14.0 mm	12.0 mm	2	PMMA	1 pc. (sterile)

### Advantages

- Circular expansion and stabilization of the capsular bag
- Reduced risk of capsular fibrosis
- Reduced risk of capsular shrinking
- Safe IOL centration in eyes with zonular dehiscence
- Stabilized conditions during phaco surgery
- CE certified

## Indications

- Anticipated zonular dehiscence
- Lens subluxation
- Expected complications during phaco surgery
- Implantation of foldable IOLs
- High risk myopic eyes

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)