



## Fluoro Touch (Fluorescein Sodium Ophthalmic Strips)

Most commonly used diagnostic strip to visualize defects or abrasions in the corneal epithelium by staining the areas of cellular loss.



### Product Description

Product Code	Description	Qty. per box
MIPL/A4	Fluorescein Sodium Ophthalmic Strips	100 strips (Sterile)
MIPL/A4	Fluorescein Sodium Ophthalmic Strips	300 strips (Sterile)

### Indications

- Useful for evaluating hard contact lens fitting
- Applanation tonometry
- Evaluating tear film breakup time
- To visualise defects in the corneal epithelium

## Specifications

- Available as sterile strips in individual pouches
- Each strip is impregnated with 1.0 mg of Fluorescein sodium
- Innovative peel open packing with transparent film on one side
- CE certified variant also available for evaluating contact lens fitting

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)