


## Micro Incision Blade

Nanoedge micro-incision phaco blades come in a wide range of sizes and shapes. These are designed to make an incision on the cornea and sclera during an ophthalmic surgical procedure.



## Lance Tip Blade

These blades are designed for making corneal/limbal incision. They are found to be very suitable for making side port.



Product Code	Description	Qty. per box
MIPL/C3/LANCE/15°	15 Degree	10 Pcs. (Sterile)
MIPL/C3/LANCE/22.5°	22.50 Degree	10 Pcs. (Sterile)
MIPL/C3/LANCE/30°	30 Degree	10 Pcs. (Sterile)
MIPL/C3/LANCE/45°	45 Degree	10 Pcs. (Sterile)

## Slit Blade


Precise width keratomes are designed to give incision of required width for insertion of phaco tip, irrigation/aspiration tip.



Product Code	Description	Qty. per box
MIPL/C3/SLIT/1.2-A	1.2 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/1.3-A	1.3 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/1.4-A	1.4 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/1.5-A	1.5 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/1.8-A	1.8 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.0-A	2.0 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.2-A	2.2 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.3-A	2.3 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.4-A	2.4 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.5-A	2.5 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.65-A	2.65 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.75-A	2.75 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/2.80-A	2.80 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/3.00-A	3.00 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/3.20-A	3.20 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/3.50-A	3.50 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/3.80-A	3.80 mm	10 pcs. (Sterile)
MIPL/C3/SLIT/4.00-A	4.00 mm	10 pcs. (Sterile)

## MVR Blade

These blades are used for making side port in small incision cataract surgery and for making precise water tight port in vitreoretinal surgery.



Product Code	Description	Qty. per box
MIPL/C3/MVR/19G	19 G Straight	10 Pcs. (Sterile)
MIPL/C3/MVR/20G	20 G Straight	10 Pcs. (Sterile)
MIPL/C3/MVR/23G	23 G Straight	10 Pcs. (Sterile)
MIPL/C3/MVR/25G	25 G Straight	10 Pcs. (Sterile)

Product Code	Description	Qty. per box
MIPL/C3/MVR/19G-A	19 G Angled	10 Pcs. (Sterile)
MIPL/C3/MVR/20G-A	20 G Angled	10 Pcs. (Sterile)
MIPL/C3/MVR/23G-A	23G Angled	10 Pcs. (Sterile)
MIPL/C3/MVR/25G-A	25G Angled	10 Pcs. (Sterile)

## Crescent Blade

These blades are designed for making scleral tunnel (pocket or groove).



Product Code	Description	Qty. per box
MIPL/C3/CRESCENT/2.0	2.0 mm Straight	10 Pcs. (Sterile)
MIPL/C3/CRESCENT/2.0-A	2.0 mm Angled	10 Pcs. (Sterile)
MIPL/C3/CRESCENT/2.2	2.2 mm Straight	10 Pcs. (Sterile)
MIPL/C3/CRESCENT/2.2-A	2.2 mm Angled	10 Pcs. (Sterile)
MIPL/C3/CRESCENT/2.5	2.5 mm Straight	10 Pcs. (Sterile)
MIPL/C3/CRESCENT/2.5-A	2.5 mm Angled	10 Pcs. (Sterile)

## Enlarger Blade

These blades are designed to give incision of required width for IOL Enlarger.



Product Code	Description	Qty. per box
MIPL/C3/ENLARGER/5.00-A	5.00 mm Angled	10 Pcs. (Sterile)
MIPL/C3/ENLARGER/5.20-A	5.20 mm Angled	10 Pcs. (Sterile)
MIPL/C3/ENLARGER/5.50-A	5.50 mm Angled	10 Pcs. (Sterile)

### Features

- Satin finish blades for minimum reflection
- Extremely sharp cutting edge
- Polycarbonate handle
- Clear polycarbonate tray
- CE certified

### Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

**Sterilization** Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com  
Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)

Ref. No.: MIPL/Pub/SPB/C3/01, Issue No.: 01, Issue Date: 01.04.2020, Rev. No. 02, Rev. Date: 01.01.2023