

## Temporary Keratoprosthesis

This disposable device is used in a vitreo retinal surgery with opaque or damaged cornea.



### Product Description

Product Code	Description	Qty. per box
MIPL/D6/7.5	Temporary Keratoprosthesis (7.5mm)	1 Pc. (Sterile)
MIPL/D6/9.0	Temporary Keratoprosthesis (9.0mm)	1 Pc. (Sterile)

### Features

- Six suture holes around the rim for creating a sealed system for the VR procedure
- The device acts as a temporary cornea during the surgery
- Biconcave design gives a widefield view of the retina
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)