

YAG Capsulotomy Lens

This lens is used for YAG laser capsulotomy procedures.



Product Description

Product Code	Description	Qty. per box
MIPL/J1	YAG Capsulotomy Lens	1 Pc.

Features

- Used for ND-YAG laser posterior capsulotomy and anterior vitreolysis.
- The +30 D lens condenses the laser energy.
- The anodized light metal rim gives a firm grip on the lens.
- Made in glass.
- CE Certified.

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)