



Steady Ring (Self Stabilizing Ring)

A disposable vitrectomy ring to stabilize the vitrectomy lenses without sutures during a VR surgery.



Product Description

Product Code	Description	Qty. per box
MIPL/D13	Steady Ring	1 Pc. (Sterile)

Features

- Special 'tripod design' for better stability during vitreo retinal surgical procedures
- Eliminates the need for a skilled assistant
- No need to suture the lens to the sclera during the surgery
- Under surface of the 3 legs are ribbed for a better grip

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase-1, New Delhi -110020, India

T: (+91-11) 4270 1030, 31 | E: Sales@madhuinstruments.com | W: www.madhuinstruments.com

Customer Care No. +91 9810866176